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Mettler-Toledo  
http://www.mt.com

METTLER TOLEDO INSTRUMENTS (SHANGHAI) CO., LTD.  
Printed by MTCS Marcom 2009/01 12107264B

## General information

Sulphide electrode is suitable for measurements of sulphide concentrations in aqueous and nonaqueous systems. Sulphide electrode may also be regarded as silver electrode; it can be the indicator electrode in potential titration for  $\text{Cl}^-$ ,  $\text{Br}^-$ ,  $\text{I}^-$ ,  $\text{CN}^-$ ,  $\text{SCN}^-$ , which can react with silver ion to form precipitation.

## Specifications

- Type of ion selective membrane:  
Solid state membrane
- Membrane impedance:  $< 5 \text{ M}\Omega$
- Measurement range:  
 $10^{-6} \text{ mol/l}$  up to approx.  $1 \text{ mol/l}$  sulphide
- Measurement temperature range:  $0\sim 50^\circ \text{ C}$
- Optimum pH range: 12 to 14
- Response time:  
when changing from lower to higher concentrations below  $10^{-4} \text{ mol/l S}^{2-}$  less than 30 seconds; when changing from higher to lower concentration several minutes.
- Interfering ion:  
Anionic surfactants (detergents) must be absent. The sulphide electrode is also sensitive to the  $\text{Ag}^+$  ion. The ratio ( $\text{Ag}^+/\text{S}^{2-}$ ) should thus not be larger than 1.

## Preparatory operations

### a. Activation

Before the sulphide electrode was used, it should be soaked in  $0.01 \text{ mol/l}$  sulphide solution for half an hour and rinsed repeatedly with deionized water till blank potential is steady. The sulphide electrode should be dabbed with paper tissue.

### b. Preparation of standard solution

Weight out  $250 \text{ g Na}_2\text{S}\cdot 9\text{H}_2\text{O}$  for a  $1 \text{ mol/l S}^{2-}$  standard sulphide stock solution and dissolve in diluent in a  $1000 \text{ ml}$  calibrated flask and make up to the mark. Calibrate the stock solution by iodimetry and store it in a refrigerator. Prepare solutions with different sulphide

concentrations by diluting the standard stock solution with diluent as appropriate.

### c. Reference electrode

Remove cap from electrolyte filling port. Slightly loosen the ground-joint diaphragm (by careful turning) so that outer reference electrolyte flows out. Tighten up ground joint again (but not too tight). Top up with  $1 \text{ mol/l NaOH}$  electrolyte to filling port for outer-mediator electrolyte and recover the cap. Clutch the electrode and shake up and down in order to reduce air bubbles.

## Calibration and measurement

If you are using **ionic meter**, proceed as follows to measure sulphide ion concentration:

### a. 2 point calibration

Prepare two standard solutions by diluting the standard stock solution with diluent as appropriate. One standard has 10 times difference in concentration from another one. The standard solutions can be prepared with any concentration unit.

According to the menu setting process of your ionic meter, set the lower concentration standard as the first calibration point, and the higher one as the second calibration point.

Place the sulphide electrode and reference electrode together in the lower standard solution, and start the first point calibration.

Place two electrodes in the higher standard solution, and start the second point calibration. Ionic meter will determine the slope automatically.

### b. Measurement Sample

Rinse two electrodes with deionized water.

To  $15 \text{ ml}$  sample solution add  $10 \text{ ml SAOB}$  solution and stir. Place two electrodes in the mixed solution and start measurement.

If you are using **pH/mV meter**, proceed as follows to measure sulphide ion concentration:

### a. Plotting calibration diagram

According to pH/mV meter instruction manual, select mV mode. Place the sulphide electrode and

reference electrode in a series of sulphide solutions one after another and measure their mV from lower to higher concentrations.

Plot measured mV values on a calibration diagram against the logarithm of the sulphide concentration of the corresponding solution (use logarithm paper).

### b. Measurement Sample

Rinse two electrodes with deionized water after each measurement.

To  $15 \text{ ml}$  sample solutions add  $10 \text{ ml SAOB}$  solution and stir. According to above method measure the mV of sample solution. Use the calibration diagram to read off the sulphide concentration of the sample solutions from the measured mV values.

It is advisable when carrying out lower sulphide concentration determination to apply a standard addition method.

## Maintenance

Deposits on the crystal membrane can be removed by carefully wiping with a damp cellulose tissue. Take care to avoid damages (scratches) to the crystal surface. Afterwards, conditioning in dilute sulphide solution (approx.  $0.01 \text{ mol/l}$ ) during 5 to 10 minutes is necessary.

## Storage

After the sulphide electrode was used, it should be rinsed to blank potential and then wiped off. If not used for a long time, it may be stored dry in the protective tube.

## Accessories

LE302 Ref. electrode  
(Order No. 12107202 or 12107204)

Diluent:

Prepare as follows: weigh out  $40 \text{ g NaOH}$  and  $5 \text{ g ascorbic acid}$ , then dissolve in deionized water and make up to 1 liter.


Sulphide Antioxidant Buffer (SAOB):

Prepare as follows: weigh out  $10 \text{ g NaOH}$  and

$1.25 \text{ g ascorbic acid}$ , then dissolve in deionized water and make up to  $100 \text{ ml}$ .

Intermediate electrolyte for the reference electrode:  $1 \text{ M NaOH}$

**Order No.: 12107250 or 12107260**

- 1 Electrode cap
- 2 Plastic electrode body
- 3 Internal lead out
- 4 Electric screw
- 5 Ion-selective membrane
- 6 Filling port for inner electrolyte
- 7 Filling port for outer electrolyte
- 8 Internal reference electrolyte
- 9 Ag/AgCl reference element
- 10 Porous ceramic diaphragm
- 11 Outer reference electrolyte
- 12 Ground-joint diaphragm