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Mettler-Toledo
<http://www.mt.com>

METTLER TOLEDO INSTRUMENTS (SHANGHAI) CO., LTD.
Printed by MTCS Marcom 2009/01 12107132B

General information

Potassium electrode is intended for measuring the potassium content of aqueous system. Its preferred applications are to be found in the fields of food-products, soil, fertilizer, beverage and biological fluid (blood, urine, perspiration, etc.) analysis.

Specifications

- Type of ion selective membrane: PVC polymer-matrix
- Membrane impedance: < 10MΩ
- Measurement range: 10⁻⁶mol/l up to 1mol/l K⁺
- Operating temperature range: 0~50° C
- Optimum pH range: 2 to 12

- Response time: when passing from lower to higher concentrations, a few seconds.

Interfering ion: cationic surface-tension agents (detergents) must be absent; the molar proportions of the following ions to the potassium ion must lie below the following levels (ratios): Rb⁺/K⁺<5×10⁻³, Cs⁺<0.02, NH₄⁺<1, Ba²⁺, Sr²⁺<20, Ca²⁺<500, H⁺<1000, Na⁺, Mg²⁺<2000.

Preparatory operations

a. Filling internal solution

Before the potassium electrode was used, place electrode with the head downwards and unscrew module from the electrode body. Fill electrode body and membrane module completely with internal filling solution ensuring that there are no air bubbles. Screw slowly the membrane module right home and then immerse it in 0.01 mol/l potassium solution for an hour and rinsed repeatedly with deionized water till blank potential is steady.

b. Preparation of standard solution

Dry AR grade potassium chloride at 120° C for 2 hours. Weight out 7.456g KCl (for a 0.1 mol/l K⁺ solution), or 1.907g KCl (for a 1.0 g/l K⁺ solution), dissolve in deionized water in a

1000ml calibrated flask and make up to the mark. Prepare solutions with different potassium concentrations by diluting the standard stock solution with deionized water as appropriate.

c. Reference electrode

Remove cap from electrolyte filling port. Slightly loosen the ground-joint diaphragm (by careful turning) so that out reference electrolyte flows out. Tighten up ground joint again (but not too tight). Top up with 0.9 mol/l Al₂(SO₄)₃ electrolyte to filling port for outer electrolyte and recover the cap. Clutch the electrode and shake up and down in order to reduce air bubbles.

Calibration and measurement

If you are using **ionic meter**, proceed as follows to measure potassium ion concentration:

a. 2 point calibration

Prepare two standards; one standard has 10 times difference in concentration from another one. The standard solutions can be prepared with any concentration unit.

According to the menu setting process of your ionic meter, set the lower concentration standard as the first calibration point, and the higher one as the second calibration point.

To 100ml the lower concentration standard add 2ml ISA solution, and stir. Place the potassium electrode and reference electrode together in the mixed solution, start the first point calibration.

To 100ml the higher concentration standard add 2ml ISA solution, and stir. Place two electrodes in the mixed solution, start the second point calibration.

After 2 points calibration, the meter automatically determines the calibration slope.

b. Measurement Sample

Rinse two electrodes with deionized water.

To 100ml sample solution add 2ml ISA solution and stir. Place two electrodes in the mixed solution and start measurement.

If you are using **pH/mV meter**, proceed as follows to measure potassium ion concentration:

a. Plotting calibration diagram

Select two or more standard solutions including the sample's potassium concentration. The standard solutions can be prepared with any concentration unit. According to pH/mV meter instruction manual, select mV mode.

To 100ml standard solutions add 2ml ISA solution and stir, respectively. Place the potassium electrode and reference electrode in the mixed solution and measure the mV of the mixed solution from lower to higher concentrations.

Plot measured mV values on a calibration diagram against the logarithm of the potassium concentration of the corresponding solution (use loggraph paper).

b. Measurement Sample

Rinse two electrodes with deionized water after each measurement.

To 100ml sample solutions add 2ml ISA solution and stir. According to above method measure the mV of sample solution. Use the calibration diagram to read off the potassium concentration of the sample solutions from the measured mV values.

It is advisable when carrying out lower potassium determination to apply a standard addition method.

Maintenance

Deposits on the active surface of the membrane module may be removed by immersing the potassium electrode in deionized water or in a diluted solution of an acid or a base (<0.01 mol/l) for a few minutes. After rinsing, recondition the electrode by immersion for a few hours in 0.01 mol/l standard K⁺ solution.


Storage

After measurement, potassium electrode must be rinsed. The potassium electrode should be stored dry in the air, preferably in the protective tube supplied with the electrode.

Accessories

LE302 Ref. electrode
(Order No. 12107202 or 12107204)
Bridge electrolyte: 0.9mol/l Al₂(SO₄)₃
Ionic Strength Adjuster (ISA): 0.9mol/l Al₂(SO₄)₃

Order NO: 12107130 or 12107100

- 1 Electrode cap
- 2 Plastic electrode body
- 3 Internal lead out
- 4 Electric screw
- 5 Ion-selective membrane
- 6 Filling port for inner electrolyte
- 7 Filling port for outer electrolyte
- 8 Internal reference electrolyte
- 9 Ag/AgCl reference element
- 10 Porous ceramic diaphragm
- 11 Outer reference electrolyte
- 12 Ground-joint diaphragm